

Fri Sep 08 2023 10:37:50 GMT-5 OPENSIGNALS VERSION: Public Build 2022-05-16

Author: PLUX

# **Devices**

98:D3:51:FD:71:01			address 98:D3:5	51:FD:71:01	type bitalino_rev	sampling rate 100 Hz
channel label	sensor	resolution	channel	label	sensor	resolution
A1 A1	EMGBITREV	10	П	INPUT	BTN	1
			12	INPUT	Generic Input	1
			01	OUTPUT	LED	1
			02	OUTPUT	BUZ	1
			DAC	DAC	PWM	8



00:50.0

00:60



00:30.0

00:40.0

# Processed channels

00:10.0

### 98:D3:51:FD:71:01

A1

0.00:00

FILTER DETECTION OPTIONS

00:20.0

applied bandpass - [10,100] Hz cutoff - Order: 6 applied to TKE: smooth window of 0.5 seconds

# Muscle Activation results

## **Activations**

find all muscle activations

Diff	-	_							-	-			
Ratio	_	_							-	-			
Glossaru	J		-	-	-	-	-	-			-	-	-

1. MVC: Maximum Voluntary Contraction.

2. TKE: Teager-Kaiser Energy operator.

